





Priority Principle™: Brachial Plexus Injury



Priority	Principle	Frequency	Time
-2	<48 Hours (Acute Injury)	5-1000 Hz	3 Minutes

# of treatments					

Priority Principle™: Brachial Plexus Injury




Priority	Principle	Frequency	Time
0	Pain (PRN)	(PRN)	
2nd	Inflammation	5-250 Hz	DOSE
3rd	Spasms	1000 Hz PRN	Pontinen's Principle
4th (a)	Pain (Systemic)	1000 Hz NRT	See Method
OR			
(b)	Pain (Local)	1000, 3000 or 5000 Hz	3-5 minutes

# of treatments					




Priority Principle™: Brachial Plexus Injury



Priority	Principle	Frequency	Time
5th	(a) Tissue Repair (Primary)	5-500 Hz	DOSE
	and		
	(b) Tissue Repair (Secondary)	50 Hz PHT @ subclavian artery	5 minutes


# of treatments					

Priority Principle™: Brachial Plexus Injury



Priority	Principle	Tendinitis	Time
6th	ROM	1000, 3000 or 5000 Hz	1-2 minutes per point
7th	Functional Strength	5-250 Hz	1-2 minutes per point

Priority Principle™: Brachial Plexus Injury



Priority	Principle	Frequency	Time
-2	<48 Hours (Acute Injury)	5-1000 Hz	3 Minutes
0	Pain (PRN)	(PRN)	
2nd	Inflammation	5-250 Hz	DOSE
3rd	Spasms	1000 Hz	Pontinen's Principle
4th (a)	Pain (Systemic)	1000 Hz NRT	See Method
OR			
(b)	Pain (Local)	1000, 3000 or 5000 Hz	3-5 minutes
5th (a)	Tissue Repair (Primary)	5-500 Hz	DOSE
and			
(b)	Tissue Repair (Secondary)	50 Hz PHT (subclavian)	5 minutes
6th	ROM	1000, 3000 Hz	2-3 minutes
7th	5-250 Hz	1-2 minutes per point	

