



# Priority Principle™: Calcific tendinitis




Priority	Principle	Frequency	Time
3rd	Spasms	1000 Hz PRN	Pontinen's Principle
4th	(a) Pain (Systemic)	1000 Hz NRT	See Method
	OR		
	(b) Pain (Local)	1000, 1000-3000 or 5000 Hz	3-5 minutes

# Priority Principle™: Calcific tendinitis



Priority	Principle	Frequency	Time
6th	ROM	1000, 3000 Hz	2-3 minutes
7th	Functional Strength	1-250 Hz	1-2 minutes per muscle

# Priority Principle™: Calcific tendinitis



	Priority	Principle	Calcific tendinitis	Time	
Cryotherapy, Ultrasound, Iontophoresis, Extracorporeal shock wave therapy (ECSW)	3rd	Spasms	1000 (PRN) with Photoprobes	Pontinen's Principle	
	4th	(a)	Pain (Systemic)	1000 Hz NRT	3-5 minutes
		(b)	Pain (Local)	1000, 1000-3000 or 5000 Hz	1-2 minutes per point
	6th	ROM	1000 or 5000 Hz	1-2 minutes per point	
	7th	Functional Strength	1-250 Hz	1-2 minutes per muscle	