


Priority Principle™: Lumbar Osteoarthritis




	Method	Dose	Notes
	Chronic	50 Hz x 3-5 minutes	URP, or PHT (shown) 50 Hz PHT (abdominal aorta)
	Inflammation	1-250 Hz or 250 Hz x 2 minutes	TARGET/DOSE preferred
	Spasm	1000 Hz x 2 minutes per spasm	Pontinen's Principle
	Pain	1000-3000 Hz for 3 minutes -----or 3000 Hz x 2-5 minutes	Systemic target (NRT) or Lumbar Sympathetic Plexus Laser Block -----or Local targets
	Tissue Repair	500-1000 Hz for 2 minutes or TARGET/DOSE -----and 50 Hz x 5 minutes	Local targets -----and Systemic target (PHT) at Abdominal aorta
	ROM	1000 Hz x 1 minute	@2-3 location around the joint line

Priority Principle™: Lumbar Osteoarthritis



Method	Dose	Notes
	50 Hz x 3-5 minutes	URP, or PHT (shown) 50 Hz PHT (abdominal aorta)

Priority Principle™: Lumbar Osteoarthritis




Method	Dose	Notes
 Inflammation	Inflammation	1-250 Hz or 250 Hz x 2 minutes
 Spasm	Spasm	1000 Hz x 2 minutes per spasm PRIN, palpate for appropriate locations
 Pain	Pain	1000-3000 Hz for 3 minutes -----or 3000 Hz x 2-5 minutes

Priority Principle™: Lumbar Osteoarthritis



Method	Dose	Notes
	500-1000 Hz for 2 minutes or TARGET/DOSE -----and 50 Hz x 5 minutes	Local targets -----and Systemic target (PHT) at Abdominal aorta

Priority Principle™: Lumbar Osteoarthritis



Method	Dose	Notes
	1000 Hz x 1 minute	@2-3 location around the joint line