


Priority Principle™: Local Pain - Options and Settings



Priority	Principle	Setting	Notes
	At site of Acute pain	1000, 3000, 5000, 1000-3000 Hz	3-5 min per location, scanning slowly no < 50 cm ²
	At site of Chronic pain	1000, 3000, 5000, 1000-3000 Hz	5-10 min per location, scanning slowly no < 50 cm ²
	Muscle (Pontinen's Principle)	1000 or 3000 Hz with Photoprobes	Accessory muscle spasm and trigger points
	Tendon	1000, 3000, 5000, 1000-3000 Hz	3-5 min per location
	Ligament (Joint): at any tender or AHSHI points (typically 3-4)	1000, 3000, 5000, 1000-3000 Hz	1-2 min with LaserStim or Photoprobe attachment
	Nerve	500-1000, 1000, or 1000-3000 Hz	The nerve root should be treated on the affected side for 1-2 min + 4 minutes of scanning along the affected dermatome

Pain Level	Setting
1-3	1000 Hz
4-6	3000 Hz
7-10	5000 Hz

Priority Principle™: Local Pain




Method	Dose	Notes
 Acute pain	1000, 3000, 5000, 1000-3000 Hz	At site of Acute pain 3-5 min per location, scanning slowly no < 50 cm ²
 Chronic pain	1000, 3000, 5000, 1000-3000 Hz	At site of Chronic pain 5-10 min per location, scanning slowly no < 50 cm ² ,

Pain Level	Setting
1-3	1000 Hz
4-6	3000 Hz
7-10	5000 Hz

Priority Principle™: Local Pain - Muscular




- Supraspinatus (1)
- Rotator Cuff (2)
- Epicondylitis Medial (3)
- Epicondylitis Lateral (4)
- Extensor Pollicis Longus (5)
- Hip Flexor (6)
- Adductor (7)
- Hamstring (8)
- Trigger Finger (9)
- Patellar (10)
- Achilles' (11)
- Plantar Fascia (12)




Method	Dose	Notes
	1000, 3000, 5000, 1000-3000 Hz	3-5 min per location


Priority Principle™: Local Pain – Nerve Root and Trunks



Method	Dose	Notes
	500-1000, 1000, or 1000-3000 Hz	the nerve root should be treated on the affected side for 1-2 min + 4 minutes of scanning along the affected dermatome




Priority Principle™: Local Pain



Method	Dose	Notes
	1000, 3000, 5000, 1000-3000 Hz	1-2 min with Photoprobe attachment At any painful joints (typically 3-4 tender points or AHSIH points)

No of AHSIH Points	Affected joint
2-3	Fingers
3-4	Wrist
4	Elbow
4-5	Shoulder
3-4	Cervical Spine
3-4	Thoracic Spine
3-4	Lumbar Spine
4-5	Hip
3-4	Knee
3-4	Ankle
2-3	Toes

Priority Principle™: Local Pain - Spasm



Method	Dose	Notes
 Spasm	1000 or 3000 Hz with Photoprobes	Pontinen's Principle Spasms: at palpable spasms in affected area, active and latent

	Corporal Probe: Corporal acupuncture is done on points on the body. This probe can also be used for muscle trigger and motor points.
	Utility Probe: This is an inter cavity probe designed for treatment of the mouth and ears. It may additionally be used as a point probe for patients that may find the corporal probe uncomfortable.