





# Priority Principle™: Knee Sprain (Medial Shown)



Priority	Principle	Frequency	Time
-2	<48 Hours (Acute Injury)	5-1000 Hz	3 Minutes

# of treatments					

# Priority Principle™: Knee Sprain (Medial Shown)



Priority	Principle	Frequency	Time
0	Pain (PRN)	(PRN)	
1st	Swelling	1000, 3000, 1000-3000 1-2 min per location	Oshiro's Principle, Proximal Priority Principle
2nd	Inflammation	50 Hz	DOSE
3rd	Spasms	1000 Hz PRN	Pontinen's Principle
4th	(a) Pain (Systemic)	1000 Hz NRT	See Method
	(b) Pain (Local)	1000, 3000 or 5000 Hz	3-5 minutes

# of treatments					

# Priority Principle™: Knee Sprain (Medial Shown)




Priority	Principle	Frequency	Time
5th	(a) Tissue Repair (Primary)	5-250 Hz	DOSE
	and		
	(b) Tissue Repair (Secondary)	50 Hz PHT (Anterior through abdomen)	5 minutes

# of treatments					




# Priority Principle™: Knee Sprain (Medial Shown)



Priority	Principle	Frequency	Time
6th	ROM	1000, 3000 Hz	2-3 minutes

# of treatments					

# Priority Principle™: Knee Sprain (Medial Shown)



Priority	Principle	Frequency	Time
-2	<48 Hours (Acute Injury)	5-1000 Hz	3 Minutes
0	Pain (PRN)	(PRN)	
1st	Swelling	1000, 3000, 1000-3000 1-2 min per location	Oshiro's Principle, Proximal Priority Principle
2nd	Inflammation	50 Hz	DOSE
3rd	Spasms	1000 Hz	Pontinen's Principle
4th (a)	Pain (Systemic)	1000 Hz NRT	See Method
<b>OR</b>			
(b)	Pain (Local)	1000, 3000 or 5000 Hz	3-5 minutes
5th (a)	Tissue Repair (Primary)	5-250 Hz	DOSE
and			
(b)	Tissue Repair (Secondary)	50 Hz PHT (Popliteal)	5 minutes
6th	ROM	1000, 3000 Hz	2-3 minutes

