



# Effectiveness of rehabilitation for patients with Subacromial impingement syndrome: a systematic review



Lori A Michener, PhD, PT, ATC, SCS, Matthew K Walsworth, PT, ECS, OCS, Evie N Burnet, PT  
Journal of Hand Therapy, Volume 17, Issue 2, April–June 2004, Pages 152–164  
<http://dx.doi.org/10.1197/j.jht.2004.02.004>

The limited evidence currently available suggests that exercise and joint mobilizations are efficacious for patients with SAIS.

Laser therapy appears to be of benefit only when used in isolation, not in combination with therapeutic exercise.

Ultrasound is of no benefit, and acupuncture trials present equivocal evidence.

Ultrasound is of no benefit, and acupuncture trials present equivocal evidence. The low to mediocre methodologic quality of the included studies limit these findings for the development of useful clinical practice guidelines. Further trials are needed to investigate these rehabilitation interventions, the superiority of one intervention over another, and the long-term outcomes of rehabilitation. Moreover, it is imperative that clinical guidelines are developed to indicate those patients who are likely to respond to rehabilitation.





# Red (660 nm) and infrared (830 nm) low-level laser therapy in skeletal muscle fatigue in humans: what is better?



Patrícia de Almeida, Rodrigo Álvaro Brandão Lopes-Martins, Thiago De Marchi, Shaiane Silva Tomazoni, Regiane Albertini, João Carlos Ferrari Corrêa, Rafael Paolo Rossi, Guilherme Pinheiro Machado, Daniela Perin da Silva, Jan Magnus Bjordal, Ernesto Cesar Pinto Leal Junior

*Lasers in Medical Science* 2012 March 2012, Volume 27, Issue 2, pp 453-458,

In animal and clinical trials low-level laser therapy (LLLT) using red, infrared and mixed wavelengths has been shown to delay the development of skeletal muscle fatigue. However, the parameters employed in these studies do not allow for a general conclusion regarding the best wavelength for delaying skeletal muscle fatigue. With this perspective in mind, we compared the effects of red and infrared LLLT on skeletal muscle fatigue.

**We conclude that both red than infrared LLLT are effective in delaying the development skeletal muscle fatigue and in enhancement of skeletal muscle performance.**

male volunteers. They were treated with active red LLLT, active infrared LLLT (660 or 830nm, 10 mW, 17.8J/cm<sup>2</sup>, 10 min) or placebo LLLT at four points of the biceps brachii muscle for 3 min before exercise (voluntary isometric elbow flexion). Mean peak force was significantly greater ( $p < 0.05$ ) following red (12.14%) and infrared LLLT (14.49%) than following placebo LLLT, and the mean average force was also significantly greater ( $p < 0.05$ ) following red (13.09%) and infrared LLLT (13.24%) than following placebo LLLT. There were no significant differences in mean average force or mean peak force between red and infrared LLLT. We conclude that both red than infrared LLLT are effective in delaying the development skeletal muscle fatigue and in enhancement of skeletal muscle performance. Further studies are needed to identify the specific mechanisms through which each wavelength acts.





# Additive effects of low-level laser therapy with exercise on subacromial syndrome: a randomised, double-blind, controlled trial




Seyyed Mohammad Jalil Abrisham, Mohammad Kermani-Alghoraishi, Rahil Ghahramani, Latife Jabbari, Hossein Jomeh, Maryam Zare

Clinical Rheumatology 2011 Clinical Rheumatology October 2011, Volume 30, Issue 10, pp 1341-1346

The subacromial syndrome is the most common source of shoulder pain. The mainstays of conservative treatment are rest, physical therapy and exercise. Low-level laser therapy (LLLT) has been popularized in the treatment of various musculoskeletal disorders. The aim of this study is to evaluate the additive effects of LLLT combined with exercise therapy on pain and shoulder range of motion (ROM) in patients with subacromial syndrome. We conducted a randomised clinical study of 80 patients who presented to clinic with subacromial syndrome (rotator cuff tendinopathy). Patients were divided into two groups. In group I (n = 40), patients were given laser treatment (pulsed infrared laser) and exercise therapy for ten sessions during a 10-week period. In group II, patients were given exercise therapy for the same period. Patients were evaluated for the pain with visual analogue scale (VAS) and shoulder range of motion (ROM) in an active and passive movement of flexion, abduction and external rotation before and after treatment. In both groups, significant post-treatment improvements were achieved in all parameters ( $P = 0.00$ ). In comparison between the two groups, a significant improvement was noted in all movements in group I ( $P = 0.00$ ). Also, there was a substantial difference between the groups in VAS scores ( $P = 0.00$ ) which showed significant pain reduction in group I. This study indicates that LLLT combined exercise is more effective than exercise therapy alone in relieving pain and in improving the shoulder ROM in patients with subacromial syndrome.




# Priority Principle™: Shoulder Subacromial Impingement



| Priority | Principle                | Frequency | Time      |
|----------|--------------------------|-----------|-----------|
| -2       | <48 Hours (Acute Injury) | 5-1000 Hz | 3 Minutes |

| # of treatments |  |  |  |  |  |
|-----------------|--|--|--|--|--|
|                 |  |  |  |  |  |

# Priority Principle™: Shoulder Subacromial Impingement




| Priority  | Principle       | Frequency             | Time                 |
|-----------|-----------------|-----------------------|----------------------|
| 0         | Pain (PRN)      | (PRN)                 |                      |
| 2nd       | Inflammation    | 50 Hz                 | DOSE                 |
| 3rd       | Spasms          | 1000 Hz PRN           | Pontinen's Principle |
| 4th (a)   | Pain (Systemic) | 1000 Hz NRT           | See Method           |
| <b>OR</b> |                 |                       |                      |
| (b)       | Pain (Local)    | 1000, 3000 or 5000 Hz | 3-5 minutes          |

| # of treatments |  |  |  |  |  |
|-----------------|--|--|--|--|--|
|                 |  |  |  |  |  |




# Priority Principle™: Shoulder Subacromial Impingement



| Priority | Principle                     | Frequency                     | Time      |
|----------|-------------------------------|-------------------------------|-----------|
| 5th      | (a) Tissue Repair (Primary)   | 5-250 Hz                      | DOSE      |
|          | and                           |                               |           |
|          | (b) Tissue Repair (Secondary) | 50 Hz PHT @ subclavian artery | 5 minutes |

| # of treatments |  |  |  |  |  |
|-----------------|--|--|--|--|--|
|                 |  |  |  |  |  |


# Priority Principle™: Shoulder Subacromial Impingement



| Priority | Principle                                                                         | Frequency                                           | Time                               |
|----------|-----------------------------------------------------------------------------------|-----------------------------------------------------|------------------------------------|
| 7th      | Functional Strength: of affected muscle(s) <u>during</u> training/activity period | 5-250 Hz pre activity or 1000 Hz immediately before | 2-3 min scanning with LaserShower* |

| # of treatments |  |  |  |  |  |
|-----------------|--|--|--|--|--|
|                 |  |  |  |  |  |

# Priority Principle™: Shoulder Subacromial Impingement



| Priority | Principle                                                                         | Frequency                                           | Time                               |
|----------|-----------------------------------------------------------------------------------|-----------------------------------------------------|------------------------------------|
| -2       | <48 Hours (Acute Injury)                                                          | 5-1000 Hz                                           | 3 Minutes                          |
| 0        | Pain (PRN)                                                                        | (PRN)                                               |                                    |
| 2nd ***  | Inflammation                                                                      | 50 Hz                                               | DOSE                               |
| 3rd      | Spasms                                                                            | 1000 Hz                                             | Pontinen's Principle               |
| 4th (a)  | Pain (Systemic)                                                                   | 1000 Hz NRT                                         | See Method                         |
| OR       |                                                                                   |                                                     |                                    |
| (b)      | Pain (Local)                                                                      | 1000, 3000 or 5000 Hz                               | 3-5 minutes                        |
| 5th (a)  | Tissue Repair (Primary)                                                           | 5-250 Hz                                            | DOSE                               |
| and      |                                                                                   |                                                     |                                    |
| (b)      | Tissue Repair (Secondary)                                                         | 50 Hz PHT (subclavian)                              | 5 minutes                          |
| 7th ***  | Functional Strength: of affected muscle(s) <u>during</u> training/activity period | 5-250 Hz pre activity or 1000 Hz immediately before | 2-3 min scanning with LaserShower* |

