




Priority Principle™: Piriformis Syndrome



Priority	Principle	Frequency	Time
0	Pain (PRN)	(PRN)	
2nd	Inflammation	50 or 5-250 Hz	DOSE
3rd	Spasms	1000 Hz PRN	Pontinen's Principle
4th (a)	Pain (Systemic)	Sacral Sympathetic Plexus Laser Block	See Method (treat affected side)
OR			
(b)	Pain (Local)	1000, 3000 or 5000 Hz	3-5 minutes

# of treatments					




Priority Principle™: Piriformis Syndrome



Priority	Principle	Frequency	Time
5th	(a) Tissue Repair (Primary)	250, 5-250 or 5-500 Hz	DOSE
	and		
	(b) Tissue Repair (Secondary)	50 Hz PHT (Anterior through abdomen)	5 minutes

# of treatments					

Priority Principle™: Piriformis Syndrome



Priority	Principle	Frequency	Time
0	Pain (PRN)	(PRN)	
2nd	Inflammation	50 or 5-250 Hz	DOSE
3rd	Spasms	1000 Hz	Pontinen's Principle
4th	(a) Pain (Systemic)	Sacral Sympathetic Plexus Laser Block	See Method (treat affected side)
	OR		
	(b) Pain (Local)	1000, 3000 or 5000 Hz	3-5 minutes
5th	(a) Tissue Repair (Primary)	250, 5-250 or 5-500 Hz	DOSE
	and		
	(b) Tissue Repair (Secondary)	50 Hz PHT (Anterior through abdomen)	5 minutes

