



Priority Principle™: Lumbar Osteoarthritis (bilateral)




Priority	Principle	Frequency	Time
-1	General Stimulation (Systemic Conditions)	URP, or PHT (shown) 50 Hz PHT (Anterior through abdomen)	3 min each

# of treatments					



Priority Principle™: Lumbar Osteoarthritis (bilateral)




Priority	Principle	Frequency	Time
0	Pain (PRN)	(PRN)	
2nd	Inflammation	5-250 or 250 Hz	DOSE
4th (a)	Pain (Systemic)	Lumbar Sympathetic Plexus Laser Block	See Method
OR			
(b)	Pain (Local)	1000, 3000 or 5000 Hz	3-5 minutes


# of treatments					



Priority Principle™: Lumbar Osteoarthritis (bilateral)




Priority	Principle	Frequency	Time
5th	(a) Tissue Repair (Primary)	250, 5-250 or 5-500 Hz	DOSE
	and		
	(b) Tissue Repair (Secondary)	50 Hz PHT (Anterior through abdomen)	5 minutes



# of treatments					




Priority Principle™: Lumbar Osteoarthritis (bilateral)



Priority	Principle	Frequency	Time
6th	ROM	1000, 3000 Hz	2-3 minutes

# of treatments					

Priority Principle™: Lumbar Osteoarthritis (bilateral)



Priority	Principle	Frequency	Time
-1	General Stimulation (Systemic Conditions)	URP, or PHT (shown)	See method
0	Pain (PRN)	(PRN)	
2nd	Inflammation	5-250 or 250 Hz	DOSE
4th (a)	Pain (Systemic)	Lumbar Sympathetic Plexus Laser Block	See Method
OR			
(b)	Pain (Local)	1000, 3000 or 5000 Hz	3-5 minutes
5th (a)	Tissue Repair (Primary)	500-1000 Hz	DOSE
and			
(b)	Tissue Repair (Secondary)	50 Hz PHT	5 minutes
6th	ROM	1000, 3000 Hz	2-3 minutes