



Priority Principle™: Spasms/Muscular Pain/Torticollis

| | | | |
|---|--|----------------------------------|----------------------|
| 3 | Spasms: at palpable spasms in affected area, active and latent | 1000 or 3000 Hz with Photoprobes | Pontinen's Principle |
|---|--|----------------------------------|----------------------|

| | |
|---|--|
|  | <p>Corporal Probe: Corporal acupuncture is done on points on the body. This probe can also be used for muscle trigger and motor points.</p> |
|  | <p>Utility Probe: This is an inter cavity probe designed for treatment of the mouth and ears. It may additionally be used as a point probe for patients that may find the corporal probe uncomfortable.</p> |

